




ObaCell[®] Fat-On-A-Chip Kit

Catalog # OS-305-01 and OS-306-01



Created with Biorender.com

Protocol Ver 2 February 2024

For Research Use Only. Not for Use in Diagnostic Procedures.


Table of Contents

| | |
|--|-----------|
| Protocol Summary | 3 |
| ObaCell® Fat-On-A-Chip Culture Setup and Maintenance | 3 |
| ObaCell® Fat-On-A-Chip Culture Harvest and Endpoints | 4 |
| General Information | 4 |
| Supplied Materials | 5 |
| Recommended Materials (Not included in the kit) | 5 |
| ObaCell® Fat-On-A-Chip Culture Setup (96-well plate) | 6 |
| ObaCell® Fat-On-A-Chip Maintenance | 7 |
| ObaCell® Fat-On-A-Chip Endpoint Assays | 8 |
| Appendix A: ObaCell® Troubleshooting | 9 |
| Appendix B: Suggested Plate Layout for the ObaCell® Fat-On-A-Chip Kit | 11 |
| Appendix C: General Guidelines and Tips for Success | 11 |
| General Guidelines | 11 |
| Tips for Success | 11 |


Protocol Summary

ObaCell® Fat-On-A-Chip Culture Setup and Maintenance


Step 1-Initial Handling




Step 2-Seeding



Step 3-Maintenance



ObaCell® Fat-On-A-Chip Culture Harvest and Endpoints



General Information

- Please refer to the safety data sheets provided with the kit for precautionary information.
- For a brief overview of the experimental procedure, please refer to the *ObaCell® Fat-On-A-Chip Quick Guide*.

Supplied Materials

| Catalog No | Product Name | Product Description | Storage Conditions |
|------------|--|--|---|
| OS-301 | ObaGel® | Human-derived, 2D and 3D applications | <48 h at 4°C, >3 months at -20°C, 1 yr at -80°C |
| OS-302 | ObaVate™ | ObaGel® activating agent | -20°C until prepared to use |
| OS-303 | ObaZolve™ | Dissolution reagent for cell recovery | 4°C for 8 weeks |
| OS-101 | Human Adipose-Derived Stromal/Stem Cells | Primary cells isolated from human adipose tissue | Liquid Nitrogen |
| OS-107 | Human Stromal Vascular Fraction | Primary cells isolated from human adipose tissue | Liquid Nitrogen |
| OS-310 | AdipoQual3D™ | Adipogenic differentiation medium | 4°C for 30 days after receipt |
| OS-309 | StromaQual3D™ | Pre-differentiation maintenance media | 4°C for 30 days after receipt |

Recommended Materials (Not included in the kit)

| Catalog No | Product Name | Product Description | Storage Conditions |
|------------|-----------------------------|--|--------------------|
| OS-001 | StromaQual™ | Complete stromal pre-differentiation medium | 4°C for 8 weeks |
| OS-008 | Cryopreservation Medium | Freezing solution for human ASCs | 4°C for 3 weeks |
| OS-304 | ObaFlow™ | Sample preparation solution for flow cytometry | 4°C for 8 weeks |
| OS-011 | Trypan Blue Exclusion Stain | Nuclear exclusion stain | 4°C for 1 year |
| | Phosphate Buffered Saline | 1X, Sterile | 4°C for 8 weeks |
| N/A | 96-well plate | Sterile, tissue culture plastic | |
| N/A | Tubes | Sterile 50 mL and 1.5 mL | |
| N/A | Micropipettes and tips | 1000, 100, and 10 µL | |
| N/A | Wet ice | For prolonged handling | |

ObaCell® Fat-On-A-Chip Culture Setup (96-well plate)

1. Thaw ObaGel® and ObaVate™ overnight at 4°C until they are completely thawed. Do not warm to 37°C. Keep reagents cold while in use.
2. Warm growth medium of choice to 37°C.
3. Prepare a 50 mL conical tube with 4 mL of growth medium.
4. Thaw the cryovial of cells until the moment the ice crystals disappear from the vial.
5. Transfer the cryovial to a BSL2 BSC and add the contents of the cryovial dropwise to the 4 mL of growth medium in the conical tube.
6. Homogenize the cell suspension by pipetting several times.
7. Perform a cell count according to your laboratory's standard operating procedures.
8. Centrifuge the conical tube at 1200 RPM for 5 minutes at room temperature.
9. Aspirate the supernatant from the cell pellet and resuspend the cells in ObaVate™ at 1 million cells per mL.
10. Calculate the volume of cell suspension needed. To seed the interior wells of a 96-well plate, we recommend a concentration of 5,000 cells per well (ASCs) or 50,000 cells per well (SVF) for 60 wells.
11. Transfer the 300 µL of cell suspension to 4.2 mL of ObaVate™ in a new 50 mL conical tube.
12. Add 1.5 mL ObaGel® to the ObaVate™ cell suspension. The total volume of the ObaGel®: ObaVate™ cell suspension should be 6 mL for a 1:3 ObaGel®: ObaVate™ mixture.
13. Pipette the solution until it is homogenous with a p1000 pipette.
14. Pipette 100 µL of the solution into the inner 60 wells of a 96-well plate (not provided) with a p200 micropipette, pipetting the mix between replicates to maintain homogeneous suspension.¹
15. Add 100 µL of sterile 1X PBS (not provided) to the outer 36 wells of the 96-well plate to mediate the effects of evaporation on long-term cultures.
16. Transfer the plate to a humidified 5% CO₂, 37°C incubator. Observe the plate after 1 hour for initial signs of gelation and again after 24 hours to observe cell growth.

Table 1. ObaCell® culture critical values.

| Cell Type | Total cells needed* | Volume of cell suspension needed | Cells per well (in ObaCell® culture) | Volume of ObaGel®:ObaVate™ cell suspension per well (1:3)** | Total number of wells seeded |
|-----------|---------------------|----------------------------------|--------------------------------------|---|------------------------------|
| ASCs | 300,000 | 300 µL | 5,000 | 100 µL | 60 |
| SVF | 3,000,000 | 300 µL | 50,000 | 100 µL | 60 |

*For 60 wells, the total number of cells needed is 300,000 (ASCs) or 3,000,000 (SVF)

**Obatala recommends accounting for pipetting error by increasing the volume by 10%.

¹ Keep suspension on ice while in use.

ObaCell® Fat-On-A-Chip Maintenance

1. Place StromaQual3D™ and AdipoQual3D™ on a benchtop to equilibrate to room temperature.
2. Examine the ObaCell® culture under a brightfield microscope for confluency, network formation, and stability.
3. Transfer the ObaCell® plate and reagents to a BSL2 BSC.
4. A plate layout is provided in Appendix B. As described in the plate layout, add 50 µL of StromaQual3D™ to three of the wells and 50 µL of AdipoQual3D™ to the other 57 wells. It is critical to add the media slowly and to the corner of each well. Adding media directly to the constructs could compromise their integrity.
5. Add additional sterile 1X PBS to the outer wells if needed.
6. Return the plate to a humidified 5% CO₂, 37°C incubator.
7. For the following media change, remove 30 µL of the spent media from each of the wells by tilting the plate and removing from the corner of the well.
8. Add 50 µL of fresh StromaQual3D or AdipoQual3D to their designated wells. It is critical to add the media slowly and to the corner of each well. Adding media directly to the constructs could compromise their integrity.
9. Steps (7) and (8) will occur 2x per week for 2 weeks. We recommended extending the culture time to 3 weeks of differentiation if you are executing metabolic assays such as lipolysis and glucose uptake to ensure cultures attain a homogeneous adipocyte population.

Table 2. Sample ObaCell® culture schedule.

| Day | Week 1 | Week 2 | Week 3 | Week 4 (with metabolic assays) |
|-----------|-------------------------------|-------------------------------|-------------------------------|--------------------------------|
| Monday | ObaCell® Setup and Incubation | Add fresh media | Remove media, add fresh media | Remove media, add fresh media |
| Tuesday | | | | |
| Wednesday | | | | |
| Thursday | | | | |
| Friday | | Remove media, add fresh media | Remove media, add fresh media | Remove media, add fresh media |

Note: Extension of cultures to week 4 is optional and recommended if you are executing metabolic assays.

ObaCell® Fat-On-A-Chip Endpoint Assays

1. Place ObaZolve™ in a water bath or equivalent piece of equipment to warm to 37°C prior to use.
2. Transfer the ObaCell® plate and reagents to a BSL2 BSC.
3. Carefully remove the spent media by tilting the plate and removing ~30- 40 µL of media from each well.
4. Add 100 µL of ObaZolve™ to each well of the ObaCell® plate.
5. Return the plate to a humidified 5% CO₂, 37°C incubator and incubate for 1 hour.
6. Transfer the plate back to the BSL2 BSC and pipette each well to observe if the gel is dissolved. If resistance is observed when pipetting, incubate the plate for an additional 30 minutes.
7. Collect each well in a tube².
8. Centrifuge each tube at 1200 RPM for 5 minutes.
9. Remove the supernatant with a pipette while also ensuring that the cell pellet remains intact.
10. After resuspension of the cell pellet, proceed the following endpoints: cell passage and expansion, cryopreservation of the cells, flow cytometry setup, or preservation and storage of the cell pellet for RNA isolation.

² If any gel fragments are observed at the base of the well, wash each well with an additional volume of 1X PBS to collect the remaining cells. Transfer the wash solution to the microcentrifuge or conical tubes.

Appendix A: ObaCell[®] Troubleshooting

| Problem | Reason | Solution |
|---|---|---|
| No gelation observed | ObaGel [®] pre-warmed at 37°C or equivalent prior to use | Do not heat ObaGel [®] prior to use, maintain on ice during use |
| | Repeated freeze-thaws | Aliquot ObaGel [®] into smaller volumes and thaw as needed |
| | Improper extended storage | Store at -20°C for up to one year from receipt date. Thaw and store at 4°C for no more than 48 hr prior to use. |
| Gelation observed during pipetting or handling | ObaGel [®] pre-warmed at 37°C or equivalent prior to use | Do not heat ObaGel [®] prior to use, maintain on ice during use |
| | ObaGel [®] not maintained at 4°C or on ice prior to use | Ensure ObaGel [®] is kept at 4°C or on wet ice prior to use |
| Contraction of 3D constructs | Physical disruption of the gel | Careful handling of the plates and manipulation of the constructs. Slowly remove and add media from edge of the well, tilt the plate to reduce shear stress at the gel-liquid interface. |
| | Inhomogeneous seeding solution | Pipette ObaGel [®] /ObaVate [™] cell suspension several times immediately prior to seeding until solution appears homogenous. Pipette solution while seeding, every 3-6 replicates. |
| Rapid spheroid or organoid formation in 3D constructs | Too much media present at the gel-liquid interface | Collect conditioned media between feedings. Feed with recommended volumes. |
| | Initial seeding density of cells is too high | Reduce seeding density. Seeding density may need to be empirically determined for each cell type and lot and may depend on several factors, such as proliferation rate. |
| | Inhomogeneous seeding solution | Pipette several times to mix well and pipette several times between replicates |

| | | |
|--|--|--|
| Constructs do not appear intact, or appear fragmented in wash solution | Constructs disrupted during pipetting | Pipette slowly and carefully when adding and removing reagents. Tilt the plate to collect solution from the gel-liquid interface. Use a wider-bore pipette tip to reduce shear stress on the constructs. |
| Gel does not appear dissolved following incubation | Incubation time not long enough | Incubate for a longer period of time, up to 30 minutes for extended cultures |
| | ObaZolve™ cell recovery solution not pre-warmed to 37°C prior to use | Warm ObaZolve™ to 37°C prior to use |
| | ObaZolve™ cell recovery solution improperly stored/handled | Store ObaZolve™ at 4°C and use within 8 weeks of receipt. Purchase fresh reagents once expired. |
| No visible cell pellet | Cell pellet may have been disrupted during wash steps | Decant supernatant by inverting microtube, on final wash use a micropipette to remove any remaining volume |

Please contact Obatala Sciences technical support at support@obatalasciences.com if you have additional questions or concerns.

Appendix B: Suggested Plate Layout for the ObaCell® Fat-On-A-Chip Kit

| | | | | | | | | | | | |
|-----|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-----|
| PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS |
| PBS | CTRL | CTRL | CTRL | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | PBS |
| PBS | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | PBS |
| PBS | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | PBS |
| PBS | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | PBS |
| PBS | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | PBS |
| PBS | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | ADIPO | PBS |
| PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS | PBS |

CTRL: 3 wells of 96-well plate maintained in StromaQual3D™

ADIPO: 57 wells of 96-well plate maintained in AdipoQual3D™

PBS: Outer wells of 96-well plate maintained in sterile 1X PBS

Appendix C: General Guidelines and Tips for Success

General Guidelines

1. All personnel should be trained and certified by the Principal Investigator regarding Universal Precautions and Handling of Bloodborne Pathogens.
2. All personnel should have access to relevant safety data sheets (SDS) pertaining to the use of materials.
3. All procedures should be conducted by investigators using the appropriate personal protective equipment.
4. Any waste materials should be decontaminated and disposed of using appropriate biohazard waste containers according to laboratory standard operating procedures.
5. Wear protective eyewear during handling of cryopreserved vials of cells.
6. Conduct all procedures described in the ObaCell® Fat-On-A-Chip protocol in a BSL2 Biosafety Cabinet.

Tips for Success

1. Do not heat ObaGel® or ObaVate™ in a water bath or equivalent piece of equipment. Thaw both products overnight at 4°C.
2. Do not heat StromalQual3D™ or AdipoQual3D™ in a water bath or equivalent piece of equipment. Both products should be placed on a benchtop until they reach room temperature.

3. If there is deviation from the protocol, ensure that the ObaGel[®]:ObaVate[™] mixture is prepared at a 1:3 ratio and no more than 5,000 cells are seeded per well of a 96-well plate.
4. Do not centrifuge or vortex ObaGel[®] to remove precipitant. The precipitant is protein and will disperse upon pipetting.
5. You do not need to add culture media after initially establishing the ObaCell[®] culture. The ObaGel[®]: ObaVate[™] mixture will support the cultures for one-week with no additional media.
6. When performing media changes, remove media by tilting the plate and removing media from the corner of the well. No more than 40 μ L of spent media should be removed at a time. We suggest only removing 30 μ L of spent media per media change.
7. The following calculations are based on the seeding of 60 wells of a 96-well plate at 5,000 cells per well. Each well has a 100 μ L volume of ObaGel[®]: ObaVate[™] cell mixture. The cell mixture is composed of a 1:3 ObaGel[®]: ObaVate[™] ratio.
 - **# Wells to be seeded** = 60
 - **Total mixture needed** = 6 mL
 - **ObaVate[™] volume** = 4.5 mL
 - **ObaGel[®] volume** = 1.5 mL
 - **Cells per well** = 5,000
 - **Desired cell concentration of mixture** = 50,000 cells per mL
 - Add **300 μ L** of 1 million per mL cell solution to **4.2 mL ObaVate[™]** and **1.5 mL ObaGel[®]**